

## Deep Freezer



### Technical Specification

<b>Temperature Range</b>	: Up to -80°C
<b>Accuracy</b>	: ± 3.0°C
<b>Uniformity</b>	: ± 5.0°C
<b>Control System</b>	: PLC or PID Based

### Features :

- Constructed and fabricated on fully automatic CNC machines. Chamber Interior is fabricated of stainless steel with dull/mirror finish & exterior is of stainless steel for GMP Model /G.I. Powder coated for standard Model.
- The inner working chamber is water & vapour tight high-grade stainless steel. All seams are arc-welded & passivity for strength & durability.
- Hermetic compressors build on cascade system with CFC free refrigeration.
- Equipment mounted on heavy castor wheels for easy movement.
- Heavy duty solid stainless steel duly electro polished wire mesh type trays are provided which are easy to clean, time & efforts saving, also adjustable for convenience.

DEEP Freezer Models	Capacity	Internal Dimension In MM			External Dimension In MM			Max. No. of Trays
		Width	Depth	Height	Width	Depth	Height	
GMP-DF100L	100 Liter	450	450	610	700	820	1550	3
GMP-DF200L	200 Liter	600	600	600	850	950	1620	3
GMP-DF325L	325 Liter	600	600	900	850	950	1900	5
GMP-DF400L	400 Liter	600	630	1060	850	980	2100	6
GMP-DF600L	600 Liter	800	800	950	1050	1100	1950	5

**NOTE :** Other Sizes and Customize models of Deep Freezer are available on request.