

Hot Air Oven



Technical Specification :

Temperature Range : Ambient+10°C to 300°C
Accuracy : ± 2.0°C
Control System : Operation through digital PID controller with inbuilt cycle timer auto cutoff.

Sr. No.	Working Chamber Size (in mm)	Capacity (in Ltr.)
1	450 x 450 x 450	90
2	450 x 450 x 600	120
3	600 x 600 x 600	200
4	600 x 600 x 900	325

- Double walled interior and exterior S.S. 304 / 316 with matt finish insulated with high grade rock wool insulation.
- S.S. Wire Mesh Trays with adjustable Support.
- Forced air circulation for uniform temperature.