



Technical Specification :

Temperature Range : Ambient +10°C to 200°C

Accuracy : $\pm 2.0^{\circ}\text{C}$

Vacuum Range : More than -700mm of Hg.

Control System : Operation through digital PID controller.

- Double walled interior and exterior S.S. 316 / 304 with matt finish insulated with high grade rock wool insulation.
- Door with glass window, silicone gasket & positive latch to maintain seal at all vacuum Levels.

Round Vacuum Oven

Sr. No.	Working Chamber Size Dia x Depth (in inch)
1	9" x 12"
2	12" x 15"
3	15" x 18"

Rectangular Vacuum Oven

Sr. No.	Working Chamber Size Dia x Depth (in inch)
1	12" x 12" x 12"
2	20" x 20" x 21"