


# Web based Software (AcquiSys)



## SOFTWARE NAME : "ACQUISYS" KESAR CONTROL SYSTEMS

### Silent Features of "ACQUISYS" Software

- Meets all 21 CFR part 11 requirements
- Ethernet communication with PLC
- Secured audit trail with date and time stamp (users, equipments, all events/alarm, SMS)
- Real time analysis & graphical analysis
- All connected equipment's utility status, process parameters & active alarms can online observed by software
- Auto & manual backup with restore data facility
- No limit for user account creation
- Forced password change on first time login
- Frequency based password expiry
- Account block with 3 wrong login attempt (unlock only by administration)
- 8 to 20 alphanumeric password length
- Online downloading for all stored data of PLC, such that master data, scanner data, mkt calculation, event/alarms, master graph & scanner graph
- GMP-compliant, tamper proof reports of measured data
- E-signature Facility for all Reports
- Report header and footer with logo can customize
- Auto time synchronize facility
- Time of auto log off in idle condition of PC (10 MINS LATER)
- Password reset server email configuration
- SOFTWARE IQ/OQ Docs provided
- Department wise equipment report facility
- Suitable for each PLC & PID based equipment manufactured by **KESAR**

